



2905 VETERANS MEMORIAL HIGHWAY - SUITE #2 800.645.5520 LANGERBIOMECHANICS.COM

The OHI Family of Brands
















Additions & Modifications Specification Guide

Langer Biomechanics has been providing healthcare professionals lower extremity solutions for more than forty years. Our broad range of custom foot orthotic shells, top covers, additions and modifications allow you to address the biomechanical needs of patients suffering from lower extremity ailments. This guide was created for you to easily determine the ideal combination of orthotic components for your patient's unique needs. These additions and modifications are the most popular, we can accommodate a range of others upon your request.

Gait Plate Extension - Modification



Not available on First Choice EXCEPT Pediatric

Specifications:

Height and distance depend on the cast

Condition(s):

To correct in toe/out toe gait

Orthotics: All rigid & semi-rigid orthotics

EXCEPT Superform

Lateral Clip - Addition / Modification



Not Available for FIRSTCHOICE

Height and distance depend on the cast

Condition(s):

Ankle sprains and instability.

Orthotics: All except fashion devices

Lateral Flange - Addition



Not available on First Choice EXCEPT Pediatric

Specifications:

Height and distance depend on the cast

Condition(s):

Ankle sprains and instability.

Orthotics: All except fashion devices

Medial Flange - Addition



Not available on First Choice EXCEPT Pediatric

Specifications:

Height and distance depend on the cast Condition(s)

Dropped navicular, pes planus, stabilize the foot

Orthotics: Can be added to all orthotics

Heel Punch (with fill) - Modification



Not Available for FIRSTCHOICE

Specifications:

11/4" standard and 1" for slim devices Slow recovery PPT® or PPT®

Condition(s):

Heel pain

Orthotics: All except Slimthotic

Deep Heel Seat - Modification



Specifications:

Condition(s):

For heel pain, add medial and lateral flanges to help control the rear foot

Orthotics: Can be added to all orthotics

Morton's Extension to Shell - Modification



Specifications: Based on the cast

Condition(s):

Severe hallux rigidus

Orthotics: All except Slimthotic

1st ray cutout - Modification



Specifications:

Based on the cast

Condition(s):

Hallux Limitus

Orthotics: Can be added to all orthotics

Kinetic Wedge - Addition



Specifications: Based on the cast

Medial, lateral or custom

Condition(s):

Functional hallux limitus

Not Available for FIRSTCHOICE Orthotics: Can be added to all orthotics

PPT® Forefoot Extension to Toes - Addition



Specifications

1/16", 1/8" PPT®

Condition(s):

Extra cushioning for comfort and shock absorption

Orthotics: Can be added to all orthotics

Plastazote® Forefoot Extension to Toes - Addition



Specifications:

1/16", 1/8" Plastazote®

Condition(s):

Calluses and/or Diabetes

Orthotics: Can be added to all orthotics

PPT®/Plastazote® Forefoot Extension to Toes - Addition



3/16" total (1/16" PPT® + 1/8" Plastazote®)

Condition(s):

Calluses and/or Lesions, IPK

Orthotics: Can be added to all orthotics

Heel Spur Pad - Addition



Specifications:

1/8" PPT®

Condition(s):

Heel pain caused by spur

Orthotics: Can be added to all orthotics

Heel Cushion - Addition



Specifications 1/8" & 1/16" PPT®

Condition(s):

Heel nain

Orthotics: Can be added to all orthotics

2-4 Met Pad - Addition



Specifications:

1/8", 3/16" thickness 1/8" distal to shell

Condition(s):

Metatarsalgia

Orthotics: Can be added to all orthotics

Met Bar - Addition



Specifications

1/8" PPT® 1/4" distal

Condition(s):

1st - 5th metatarsal pain

Orthotics: Can be added to all orthotics

Toe Crest Pad - Addition



Specifications:

1/8" - 3/16" PPT®


2nd - 5th based on cast

Condition(s):

Hammer toes / claw toes

* Not Available for FIRSTCHOICE Orthotics: Can be added to all orthotics

Neuroma Pad - Addition



Specifications:

1/8" - 3/16" PPT® 1st - 4th interspace

Condition(s): Neuromas in the interspaces

Orthotics: Can be added to all orthotics

Neuroma Plug - Addition



Specifications:

1/8" - 3/16" PPT® 1st - 4th interspace

Condition(s):

Neuromas in the interspaces

Orthotics: Can be added to all orthotics

Dancer's Pad - Addition

(may be called a Reverse Morton's at other labs)



Not Available for FIRSTCHOICE

Specifications: 1/8" PPT®

Condition(s):

1st metatarsal pain/bunion

Orthotics: Can be added to all orthotics

Scaphoid Pad - Addition



Specifications

1/8" - 1/4" PPT®

Condition(s):

Medial longitudinal arch pain / support to raise the arch softly

Orthotics: Can be added to all orthotics

Morton's Extension - Addition (may be called a soft Morton's at other labs)



Specifications 1/4" crepe

Condition(s): Hallux rigidus

Not Available for FIRSTCHOICE **Orthotics:** Can be added to all orthotics

Balance Pad - Addition



Specifications 1/4" korex

1/16"- 1/4" PPT®

Condition(s): Balance out calluses or metatarsal prominences

Not Available for FIRSTCHOICE Orthotics: Can be added to all orthotics

Toe Buttress* - Addition



Specifications Plastazote®

Condition(s): Fill to accommodate amputated toes or full

transverse metatarsal amputation Orthotics: Can be added to all orthotics

Not Available for FIRSTCHOICE **Heel Lift** - Addition



Specifications:

Crepe, not recommended greater than 1/2" we suggest blending to metatarsals with added PPT® arch fill for feathering

Condition(s):

Leg length discrepancy

Orthotics: Can not be added to any device that is standard intrinsically posted (lifts sent separately)

